Boise State University

ScholarWorks

[College of Arts and Sciences Poster](https://scholarworks.boisestate.edu/as_13) **Presentations**

[2013 Undergraduate Research and Scholarship](https://scholarworks.boisestate.edu/2013_under_conf) **Conference**

4-15-2013

Chemical Warfare Agents in Plants: Biodefensive Terpenes from Sagebrush

Belinda Stierman Department of Chemistry, Boise State University

Chemical Warfare Agents in Plants: Biodefensive Terpenes from Sagebrush

Abstract

To prevent being eaten, some plants create compounds that are toxic to herbivores, and herbivores respond by creating new ways to metabolize these toxins. For example, sagebrush produces biodefensive terpenes to deter foraging by sage grouse. A recent investigation to characterize plant secondary metabolites (PSM's) in sagebrush resulted in a correlation between terpene concentration and nutritional content in local sagebrush to the habitat selection of sage grouse. This study identified many of the terpenes in sagebrush, but not all of them. A highly volatile and elusive terpene is suspected to be an important PSM that significantly affects sage grouse foraging. The purpose of the current work is to identify the structure of this PSM, as well as other unknown PSM's in sagebrush, to better understand their role as chemical warfare agents to defray foraging. Sagebrush extracts will be analyzed and components identified using Gas Chromatography-Mass Spectrometry (GC-MS).

Keywords

sagebrush, bio-defensive, terpenes, plant secondary metabolites (PSM's), sage grouse

Disciplines

Biochemistry | Plant Sciences

Abstract

To prevent being eaten, some plants create compounds that are toxic to herbivores, and herbivores respond by creating new ways to metabolize these toxins. For example, sagebrush produces bio-defensive terpenes to deter foraging by sage grouse. A recent investigation to characterize plant secondary metabolites (PSM's) in sagebrush resulted in a correlation between terpene concentration and nutritional content in local sagebrush to the habitat selection of sage grouse. This study identified many of the terpenes in sagebrush, but not all of them. A highly volatile and elusive terpene is suspected to be an important PSM that significantly affects sage grouse foraging. The purpose of the current work is to identify the structure of this PSM, as well as other unknown PSM's in sagebrush, to better understand their role as chemical warfare agents to defray foraging. Sagebrush extracts will be analyzed and components identified using Gas Chromatography-Mass Spectrometry (GC-MS).

Plants produce primary and secondary metabolites that perform an extensive array of functions. The importance of plant secondary metabolites (PSMs) in nature is grossly understated. These compounds provide many benefits for the plant including the ability to be used as a deterrent.⁽¹⁾ Sagebrush contains relatively high concentrations of PSMs; the key class of these defensive molecules are terpenes.⁽²⁾ Terpenes are the largest class of compounds in the biological volatile organic compound (BVOC) classification.⁽³⁾ Isoprene, a C_5H_8 hydrocarbon, is the basic building block of terpenes. While these compounds do have a pleasant aroma, as seen with limonene in citrus fruits, and pinene in conifers, terpenes are also toxic. (1) In the ongoing warfare between plants and herbivores, defense mechanisms are constantly evolving.⁽⁴⁾ This is observed in the phenomenon referred to as specialization. One such example can be seen in the interaction between sagebrush and sage grouse (Figure 1). (5)

Chemical Warfare Agents: Bio-defensive Terpenes from Sagebrush Belinda Stierman, † Jennifer Sorensen Forbey§ & Owen McDougal † BOISE DESTATE [†]Chemistry and Biochemistry; [§] Biological Sciences

Procedure: (6)

- 1) GC-MS analysis of terpene standards.
- 2) GC-MS Method:
	- Start temperature: 40° C, ramp $1 - 40^{\circ}$ C-125 $^{\circ}$ C at 5 $^{\circ}$ C/min,
- 3)Compare headspace to solvent extraction methods to identify terpenes
- (Table 1; Figure 4).
- them by GC-MS.

Acknowledgements

Methods

References

Future Work

Introduction

Figure 3. GC-MS instrument used to quantitate and qualitate terpenes

Conclusion

Figure 5B: Biomass was prepared by flash freezing with liquid N_2 .

Recent work has attempted to characterize the importance of terpene concentration and nutritional content in local sagebrush on the habitat selection of sage grouse.⁽⁵⁾ Though terpenes have been well researched as a whole, this study has come across an important volatile compound produced by sagebrush that significantly effects sage grouse foraging (unknown #1). The purpose of our research is to identify the structure of unknown #1, in addition to other unidentified PSMs in sagebrush, in order to elucidate the dynamics of chemical defense mechanisms on foraging. Figure 1. Sage grouse chick.

Gas Chromatography-Mass Spectrometry (GC-MS) is used because of it's accuracy in identifying volatile compounds using headspace methods (Figures 2; 3). A representative headspace vial is shown in Figure 2. ^(2,6) Traditionally a GC-MS sample is prepared by dissolving the sample in a solvent (usually methylene chloride or hexane) and injecting the solvent into the GC-MS for analysis. This method works well for non-volatile low molecular weight compounds (<400 amu), however many terpenes are BVOCs, which requires a method to isolate both volatile and the non-volatile components. A headspace sample is prepared for volatile terpenes, in a vial containing the sample in native form (solid or volatile liquid). Volatile components from complex sample mixtures can be separated from nonvolatile sample components and isolated in the headspace or gas portion of a sample vial. A sample of the gas in the headspace is then directly injected into a GC-MS for component identification. **Figure 2**. A diagram

(1) Unisicker, S.B.; Knert, G.; Gershenzon, J. "Protective Perfumes: the Role of Vegetative Volatiles in Plant

(2) Zewenger, S.; Basu, C. "Plant Terpenoids: Applications and Future Potentials" *Biotechnology and*

Table 1. A list of standard terpenes, GC retention times for headspace vs. solvent extraction methods.

Terpene	Retention time Headspace (min)	Retention Traditio (min)
a-pinene	6.35	6.32
camphene	6.77	6.75
b-pinene	7.51	7.48
a-phellandrene	8.29	8.27
1,4-cineole	8.54	8.52
p-cymene	8.82	8.78
1,8-cineole	9.03	9.01
thujone	11.16	11.14
camphor	12.36	12.29
borneol	13.09	12.99
terpineol	13.70	NA

8.52

9.01

11.14

12.29

12.99

(3) Ormeno, E.; Goldstein, A.; Niinemets, U. "Extracting and Trapping Biogenic Volatile Organic Compounds Stored in Plant Species" *Trends in Analytical Chemistry*. 2011, 30, 7, 978-989.

(4) Wittstock, U.; Gershenzon, J. "Constitutive Plant Toxins and Their Role in Defense Against Herbivores and Pathogens" *Current Opinion in Plant Biology*. 2002, 5, 4, 300-307.

(5) Frye, G.; Connelly, J.; Musil, D.; Forbey, J. Phytochemistry Predicts Habitat Selection By An Avian Herbivore at Multiple Spatial Scales. *Ecology*. 2013 (2): 308-314.

(6) Bowman, J.; et al. "Extraction Method for the Isolation of Terpenes from Plant Tissue and Subsequent Determination by Gas Chromatography" *Microchemical Journal*. 1997, 56, 10-18.

STUDENT RESEARCH INITIATIVE PROGRAM

• Headspace was superior for loading and identifying volatile terpenes by GC-MS, but it was less

Figure 5A: First, Sagebrush is harvested, then stored in -20°C freezer until ready for use.

Figure 5C: Pulverizing to a powder by mortar and pestle.

Figure 5D: Analyzed by both headspace and solvent extraction by GC-MS.

headspace (top) and traditional (bottom) methods.