Boise State University ScholarWorks

2020 Undergraduate Research Showcase

Undergraduate Research and Scholarship Showcases

4-24-2020

Faster Liposome Preparation Using Electrodialysis

Zoe E. Hutchinson Boise State University

Gamid Abatchev Boise State University

Andrew Bogard Boise State University

Jason Ward *Boise State University*

Daniel Fologea Boise State University

Faster Liposome Preparation Using Electrodialysis

Abstract

Liposomes are self-enclosed spherical lipid bilayer membranes that present a great potential for drug delivery. Liposome preparation and loading may be achieved by several methods including extrusion and sonication, all of which are rather time consuming. To alleviate this problem, we propose using electrodialysis as a fast technique for liposome preparation. Our approach is based on the fast separation of charged molecules upon their movement when exposed to electric fields. We subjected homogeneous solutions of lipids and ionic detergents to electrodialysis. The lack of amphiphilic detergents in the mixture forced the lipids to self-assemble into liposomes, incorporating dye molecules present in the initial mixture. The obtained liposomes were characterized by Dynamic Light Scattering, microscopy, and fluorescence spectroscopy. Our results show that the procedure may also be applied for preparation of liposomes capable of evading the immune system, therefore presenting potential for use as drug carriers and in vivo applications.



BOISE STATE UNIVERSITY

Abstract Liposomes are self-enclosed spherical lipid bilayer membranes that present a great potential for drug delivery. Liposome preparation and loading may be achieved by several methods including extrusion and sonication. These methods can be time consuming, require specific equipment, and of require further steps to purify produced liposomes of loaded contents. To alleviate this problem, we propose using electrodialysis as a fast technique for liposome preparation. Our approach is based on their electrodialysis as a fast technique for liposome preparation. homogeneous solutions of lipids and ionic detergents to electrodialysis through a polycarbonate membrane. The migration of the detergent out of the mixture allows the lipids to self-assemble into liposomes, incorporating designated drug/drug simulator molecules present in the initial mixture. The obtained liposomes were characterized by Dynamic Light Scattering (DLS), microscopy, and fluorescence spectroscopy (FS). Our results show that the procedure may also be applied for preparation of liposomes capable of evading the immune system, therefore presenting potential for use as drug carriers and in vivo applications with higher efficacy than traditional liposomes.

Objectives

1. Prepare liposomes using electrodialysis (ED) techniques and ionic detergents 2. Prepare and load liposomes concurrently with drug

simulating dyes

3. Characterize the prepared liposomes (e.g size, payload incorporation) using microscopy, DLS, and FS

Materials and Methods

Materials:

Lipids: Asolectin (Aso), Cholesterol (Chol), 1,2-distearoylsn-glycero-3 phosphocholine (DSPC), 1,2-distearoyl-snglycero-3-phosphoethanolamineN-[methoxy(polyethylene glycol)-2000] (DSPE-PEG)

Compositions:

A - 5 mg Aso:2 mg Chol


B - 4.1 mg DSPC:1.9 mg Chol:1.3 mg DSPE-PEG

Buffers: 135 mM KCI + 20 HEPES, 20 mM KCI + 5 mM HEPES

Detergents Used: Cholic Acid 10% W/V (CA), Sodium dodecyl sulfate 10% W/V (SDS)

Dyes: Rhodamine (Rdn) and Acridine Orange (AO).

Detergent for Release: Detergent For Release: Triton X-100 5% v/v (TX)





Faster Liposome Preparation Using Electrodialysis Zoe Hutchinson, Gamid Abatchev, Andrew Bogard, Jason Ward Mentor: Dr. Daniel Fologea

Results and Discussion

Microscopy

The liposomes were successfully prepared and loaded with ionic detergents using electrodialysis.






Before (left panels) and after (right panels) ED treatments of (2) composition A prepared with CA and 200 µM AO (3) composition A prepared with CA and 50 µM Rdn (4) composition A prepared with SDS with 35 µM R6G, all at 40x magnification.

Dynamic Light Scattering

The prepared liposomes were analyzed by DLS (Malvern Zetasizer) to characterize their hydrodynamic size and uniformity (see the examples in **Figure 5**. These liposomes were submicron size.




Assessment of Dye Loaded into Liposomes

Next, we employed fluorescence spectroscopy (Fluoromax-HORIBA) to investigate the release of incorporated dyes. The payload release was triggered by a non-ionic detergent, Triton X (Figure 6) This exploration relied on self-quenching (i.e., a significant decrease of the dye fluorescence over a certain concentration attained by the dye loaded into liposomes).


Destabilization of the liposome membrane (for example by using a detergent) induces dye leakage into the bulk, a decrease in concentration under the self-quenching level, and a substantial increase in fluorescence.



The project described was supported by the NSF Career Award 1554166 and Boise State University McNair Scholars Program. We also would like to acknowledge support from the Department of Physics and the Biomolecular Sciences Graduate Program at Boise State University.



In this example (Fig. 6), 20 µL liposomes loaded with AO were added to 1 ml of 135/20 buffer and assessed by fluorescence spectroscopy (excitation: 488 nm, emission: 525 nm). After a stable baseline was attained, addition of 50 µL Triton X 100 5% elicited a significant increase in fluorescence, indicative of dye leakage through the detergent-permeabilized membrane. This experiment indicated successful loading of AO into liposomes, similar to microscopy.



loaded during electrodialysis.

Triton x-100 suggests membrane permeabilization and release of AO

Conclusions

- Liposomes are excellent tools for incorporating water-soluble molecules into their inner cavity.

- Liposomes may be simultaneously produced, loaded, and separated from the unloaded compounds by electrodialysis in a relatively short period of time.

- Filtration and/or centrifugation may be used to further filter out aggregates and homogenize the population size.

- DSPC-Chol formulation (Composition B) was more prone to aggregation, indicative of PEG-induced stabilization over time.